



40 Terrill Park Drive, Unit 4, Concord, NH 03301  
Ph: 603-225-0180

We are liquidating used equipment in excellent shape for a great price.

All offers considered!



Equipment sold "as is" with a limited 30-day warranty



Monaco  
\$1,895



Phoenix  
\$2,995



Sorrento  
\$2,395

**Seating Matter** chairs offer safety, pressure management, and support for the user, and caregivers. These chairs are made with optimum quality and have been trialed all over the world.

**The Monaco-** reduces sliding, discomfort and shear as it helps to maintain positioning of the pelvis. Its low seat to floor height makes it easy to stand from, making it a multi-functional comfort chair recommended for mobile patients.

**The Phoenix-** provides maximum support for the upper body with a unique head support for patients with complex postural needs and spinal problems. Designed for posture care, pressure care, early mobilization and safe patient handling, the Phoenix is suitable for the most complex level of care.

**The Sorrento-** is a simple to operate and functional chair which is designed to improve and maximize comfort as well as provide postural support and pressure management. The Sorrento tilt in space chair is clinically designed for the patient but has the caregiver in mind with ease of use designed into the chair.



40 Terrill Park Drive, Unit 4, Concord, NH 03301  
Ph: 603-225-0180



**The Rowalker400-** offers safe and secure support during walking and standing training with users with impaired leg strength and balance. The Rowalker400 is an ergonomically designed walking aid for early mobilization, specifically for cardio-thoracic and abdominal post surgical users. Provides adjustable height settings, width settings, directional wheel locks, and caster locks. Weight limit is 400 pounds.

The Rowalker400  
\$2,999



The Burke Bed Tri Flex II  
\$3,995

**The Burke Bed-** The Tri-Flex II bed is multi-purpose capable. The patient surface can be raised or lowered to facilitate patient care. The bed has Trendelenburg and Reverse Trendelenburg function in the pendant as well as “Cardio-Chair” positioning. The Tri-Flex II system can also facilitate patient transport. The two foot end casters contain steering lock ability which provides directional control during transport.